


Direct Nasopharyngeal Swab Procedure




Positive



Negative



Invalid



T = COVID-19 Test Line

C = Control Line



Consult Instructions for Use



Sona Nanotech
Dartmouth, NS B2Y 4M9 Canada
Sonanano.com

External Control Dropper Bottle Procedure

1

Positive Control




Negative Control



2



15 Min



Dispense
4 Drops

Positive



C
T

Negative



C
T

Invalid



C
T



C
T

T = COVID-19 Test Line

C = Control Line



Consult Instructions for Use



Sona Nanotech
Dartmouth, NS B2Y 4M9 Canada
Sonanano.com